

DRUG OF ABUSE






One Step Drugs of abuse Test

One Step drug of abuse Test kits are lateral flow, one-step, qualitative, competitive inhibition immunoassays. The tests are based on the principle of the highly specific immunochemical reactions of antigens and antibodies. The test systems employ unique monoclonal antibodies to selectively identify specific drug of abuse and the corresponding metabolites in urine samples with a high degree of sensitivity. The cut off concentration for specific drug of abuse and its metabolites has been developed at the concentration standard set by NIDA (the National Institute on Drug Abuse) for the qualitative detection of corresponding drug in human urine specimens.

Features

- 1. Highly sensitive and specific;
- 2. Fast and easy to use, results are read visually without the need for instrumentation;
- 3. No extra reagent is needed;
- 4. Used to screen urine for the presence of corresponding drug of abuse and its metabolites;
- 5. Single test is available in both strip and cassette styles;
- 6. 2 to 10 multi-drug test panels are available (combinations subject to customers' selection).

Test Procedure





DRUG OF ABUSE



Ordering information

CAT. NO.	Description	Cut-Off	Packing
AMP 611	Amphetamine Urine Screening Strip	1000 ng/ml	100 Tests/Box; 48 Boxes/Carton
AMP 612	Amphetamine Urine Screening Cassette	1000 ng/ml	25 Tests/Box; 40 Boxes/Carton
BAR 621	Barbiturates Urine Screening Strip	300 ng/ml	100 Tests/Box; 48 Boxes/Carton
BAR 622	Barbiturates Urine Screening Cassette	300 ng/ml	25 Tests/Box; 40 Boxes/Carton
BZO 631	Benzodiazepine Urine Screening Strip	300 ng/ml	100 Tests/Box; 48 Boxes/Carton
BZO 632	Benzodiazepine Urine Screening Cassette	300 ng/ml	25 Tests/Box; 40 Boxes/Carton
COC 641	Cocaine Urine Screening Strip	300 ng/ml	100 Tests/Box; 48 Boxes/Carton
COC 642	Cocaine Urine Screening Cassette	300 ng/ml	25 Tests/Box; 40 Boxes/Carton
MDMA 651	Ecstasy Urine Screening Strip	500 ng/ml	100 Tests/Box; 48 Boxes/Carton
MDMA 652	Ecstasy Urine Screening Cassette	500 ng/ml	25 Tests/Box; 40 Boxes/Carton
MET 661	Methamphetamine (mAMP) Urine Strip	500 ng/ml	100 Tests/Box; 48 Boxes/Carton
MET 662	Methamphetamine (mAMP) Urine Cassette	500 ng/ml	25 Tests/Box; 40 Boxes/Carton
MOR 671	Morphine Urine Screening Strip	300 ng/ml	100 Tests/Box; 48 Boxes/Carton
MOR 672	Morphine Urine Screening Cassette	300 ng/ml	25 Tests/Box; 40 Boxes/Carton
MTD 681	Methadone Urine Screening Cassette	300 ng/ml	100 Tests/Box; 48 Boxes/Carton
MTD 682	Methadone Urine Screening Strip	300 ng/ml	25 Tests/Box; 40 Boxes/Carton
OPI 691	Opiates Urine Screening Strip	2000 ng/ml	100 Tests/Box; 48 Boxes/Carton
OPI 692	Opiates Urine Screening Cassette	2000 ng/ml	25 Tests/Box; 40 Boxes/Carton
PCP 6011	Phencyclidine Urine Screening Strip	25 ng/ml	100 Tests/Box; 48 Boxes/Carton
PCP 6012	Phencyclidine Urine Screening Cassette	25 ng/ml	25 Tests/Box; 40 Boxes/Carton
TCA 6021	Tricyclic Antidepressant Test Strip	1000 ng/ml	100 Tests/Box; 48 Boxes/Carton
TCA 6022	Tricyclic Antidepressant Test Cassette	1000 ng/ml	25 Tests/Box; 40 Boxes/Carton
THC 6031	Marijuana Urine Screening Strip	50 ng/ml	100 Tests/Box; 48 Boxes/Carton
THC 6032	Marijuana Urine Screening Cassete	50 ng/ml	25 Tests/Box; 40 Boxes/Carton
BUP 6041	Buprenorphine Urine Screening Strip	10 ng/ml	25 Tests/Box; 40 Boxes/Carton
KET 6051	Ketamine Urine Screening Strip	1000 ng/ml	100 Tests/Box; 48 Boxes/Carton
KET 6052	Ketamine Urine Screening Cassette	1000 ng/ml	25 Tests/Box; 40 Boxes/Carton
DOA 602 *	Multi Drug 2P Urine Screening Panel	Ref. to single test	25 Tests/Box; 40 Boxes/Carton
DOA 603 *	Multi Drug 3P Urine Screening Panel	Ref. to single test	20 Tests/Box; 33 Boxes/Carton
DOA 604 *	Multi Drug 4P Urine Screening Panel	Ref. to single test	20 Tests/Box; 33 Boxes/Carton
DOA 605 *	Multi Drug 5P Urine Screening Panel	Ref. to single test	20 Tests/Box; 33 Boxes/Carton
DOA 606 *	Multi Drug 6P Urine Screening Panel	Ref. to single test	20 Tests/Box; 33 Boxes/Carton
DOA 610 *	Multi Drug 10P Urine Screening Panel	Ref. to single test	20 Tests/Box; 33 Boxes/Carton
DOA710-715	Multi-Drug Integrated Cup (from 5-10 P)	Ref. to single test	Special packing: 100pcs/carton Size of carton: 62*42*38cm Gross weight: 11.8kg
	man brag megrated cap morns for y	rento single test	and the state of t

^{*} Combination subject to customer's selection;

^{**} Carton Outer Size: 52*45*36CM. Gross weight about 14-15.5 KG/Carton.