

FERTILITY



FSH Menopause Urine Test Kits

The FSH (human Follicular Stimulating Hormone) test results are of great clinical value to the diagnosis of many problems related with reproductive activities. One Step FSH Urine Test is a qualitative test kit for the determination of FSH concentration in urine specimens, and in turn helps the testee to get prepared both physically and mentally for the syndromes come along with menopause. Visual and qualitative test results can be read visually without any instrumental help.

Features

1. An easy to use FSH test. No additional instrument or reagent is needed;
2. High sensitivity :25MIU/ML and accuracy: 99.47%;
3. Designed for professional and home users;
4. Each test is enclosed in a moisture proof aluminum pouch with a Silicon Gel for long time storage.

Test Procedures

Strip



Cassette



Midstream



Ordering information

CAT. NO.	Description	Method	Packing	Specimen
FSH 141	FSH Urine Strip	GICA	100 Tests/Box; 48 Boxes/Carton	Urine
FSH 142	FSH Urine Cassette	GICA	25 Tests/Box; 40 Boxes/Carton	Urine
FSH 143	FSH Urine Midstream	GICA	1 test/pouch; 1200 PCS/Carton	Urine

* GICA: Lateral Flow Gold Immunochromatographic Assay;

** Carton Outer Size: 52*45*36CM. Gross weight about 14-15.5 KG/Carton.