

Fetal Doppler

FD-200 ultrasonic fetal doppler meets the fetus daily check and routine examination at home, clinic, community and hospital.



Features

- Easy and convenient operation
- High-fidelity, crystal clear sound
- Earphone and speaker are feasible
- High sensitivity doppler probe
- Low ultrasound dosage
- Display with color LCD
- **FHR digital voice funtion (option)**



Fetal Doppler



Technical Specification

- **Ultrasound frequency:** 2MHz
- **Ultrasound intensity:** <10mW/cm²
- **Power supply:**
2pcs 1.5V AA battery
2pcs DC Ni-Mh rechargeable battery
AC 220/110V,50/60Hz (for FD-200D and FD-200C)
- **Display:** 45mm×25mm LCD
- **FHR Measuring range:** 50~240bpm
- **FHR resolution:** 1bpm
- **FHR accuracy:** ±1bpm
- **Power consumption:** <1W
- **Dimension:** 135mm ×95mm ×35mm
- **Weight:** 500g

FD-200



Functions	FD-200A	FD-200D	FD-200C	FD-200B	FD-200G
Display	No	B/W	Colour	B/W	Colour
Backlight	No	Yes	Yes	Yes	Yes
Build-in speaker	Yes	Yes	Yes	Yes	Yes
Battery indicator	No	Yes	Yes	Yes	Yes
FHR display	No	Yes	Yes	Yes	Yes
FHR curve	No	No	Yes	No	Yes
FHR voice function	No	Option	Option	Option	Option
Various work modes	No	Yes	Yes	Yes	Yes
Alkalinity battery	Yes	No	No	Yes	Yes
Ni-Mh battery	No	Yes	Yes	No	No
Charge in device	No	Yes	Yes	No	No
Earphone	Option	Option	Option	Option	Option
Carry bag	Option	Option	Option	Option	Option

All specifications are subject to change without prior notice.