



Adhesion Microscope Slides



The continuous liquid-base in stains as well as chemicals and the frequent washing throughout the immunohistochemical staining workflow, especially the immediate washing after the antigen-retrieval heating do cause great risk of tissue section loss or cell loss. The surface of the CITOTEST® Color-PLUS™ microscope slides is coated with a permanent especially-formulated adhesion which reliably binds the tissue sections to the glass, considerably minimizing the possibility of tissue-section or cell loss.

The CITOTEST® Color-PLUS™ microscope slides are also highly recommended in routine H&E staining to secure the problematic tissue sections like brain and spinal cord specimens.

CAUTIONS

- 1. Store the products under dry and room temperature conditions.
- 2. Use the products as soon as possible when the package is open.
- 3. The color frosted layer may fall off the glass where caustic conditions are applied.



Recommended ★
Strongly recommended ★ ★

Applications

	PCB	PCL	PCI
HISTOLOGY AND PATHOLOGY			
Paraffin embedded section	★ ★	★ ★	★ ★
Frozen section	★ ★	★ ★	★ ★
Liquid-rich specimens	★	★	★ ★
AUTOMATED IHC STAINER			
Roche Ventana BenchMark	★	★	★ ★
Leica BOND-MAX	★ ★	★ ★	★ ★
Dako Autostainer Link 48	★ ★	★ ★	★ ★
HIER PROCESSES			
Acid Retrieval Solution	★ ★	★ ★	★ ★
Alkaline Retrieval Solution	★	★	★ ★
CYTOLOGY			
LBC, membrane	★ ★	★ ★	★ ★
LBC, centrifugation	★ ★	★ ★	★ ★
Microbiology	★	★	★
Hematology	★	★	★

ADHESION
MICROSCOPE SLIDES

Color-PLUS™ PCB Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCB series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End which is resistant to most common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils and mark pens.
Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS		BEVELED EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7141-02	0303-7141-02	0312-7141-02	0302-7141-02	0314-7141-02	0304-7141-02
Green	0313-7141-08	0303-7141-08	0312-7141-08	0302-7141-08	0314-7141-08	0304-7141-08
Purple	0313-7141-09	0303-7141-09	0312-7141-09	0302-7141-09	0314-7141-09	0304-7141-09
Orange	0313-7141-10	0303-7141-10	0312-7141-10	0302-7141-10	0314-7141-10	0304-7141-10
Pink	0313-7141-12	0303-7141-12	0312-7141-12	0302-7141-12	0314-7141-12	0304-7141-12
Tan	0313-7141-15	0303-7141-15	0312-7141-15	0302-7141-15	0314-7141-15	0304-7141-15
White	0313-7141-16	0303-7141-16	0312-7141-16	0302-7141-16	0314-7141-16	0304-7141-16
Yellow	0313-7141-17	0303-7141-17	0312-7141-17	0302-7141-17	0314-7141-17	0304-7141-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS		BEVELED EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7151-02	0303-7151-02	0312-7151-02	0302-7151-02	0314-7151-02	0304-7151-02
Green	0313-7151-08	0303-7151-08	0312-7151-08	0302-7151-08	0314-7151-08	0304-7151-08
Purple	0313-7151-09	0303-7151-09	0312-7151-09	0302-7151-09	0314-7151-09	0304-7151-09
Orange	0313-7151-10	0303-7151-10	0312-7151-10	0302-7151-10	0314-7151-10	0304-7151-10
Pink	0313-7151-12	0303-7151-12	0312-7151-12	0302-7151-12	0314-7151-12	0304-7151-12
Tan	0313-7151-15	0303-7151-15	0312-7151-15	0302-7151-15	0314-7151-15	0304-7151-15
White	0313-7151-16	0303-7151-16	0312-7151-16	0302-7151-16	0314-7151-16	0304-7151-16
Yellow	0313-7151-17	0303-7151-17	0312-7151-17	0302-7151-17	0314-7151-17	0304-7151-17

Color-PLUS™ PCL Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCL series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils, mark pens as well as most thermal transfer based slide printer, including Thermo Slide-Mate® (Thermo Fisher) slide printer, and Signature® (Primera) slide printers
Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7181-02	0303-7181-02	0312-7181-02	0302-7181-02
Green	0313-7181-08	0303-7181-08	0312-7181-08	0302-7181-08
Purple	0313-7181-09	0303-7181-09	0312-7181-09	0302-7181-09
Orange	0313-7181-10	0303-7181-10	0312-7181-10	0302-7181-10
Pink	0313-7181-12	0303-7181-12	0312-7181-12	0302-7181-12
Tan	0313-7181-15	0303-7181-15	0312-7181-15	0302-7181-15
White	0313-7181-16	0303-7181-16	0312-7181-16	0302-7181-16
Yellow	0313-7181-17	0303-7181-17	0312-7181-17	0302-7181-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7191-02	0303-7191-02	0312-7191-02	0302-7191-02
Green	0313-7191-08	0303-7191-08	0312-7191-08	0302-7191-08
Purple	0313-7191-09	0303-7191-09	0312-7191-09	0302-7191-09
Orange	0313-7191-10	0303-7191-10	0312-7191-10	0302-7191-10
Pink	0313-7191-12	0303-7191-12	0312-7191-12	0302-7191-12
Tan	0313-7191-15	0303-7191-15	0312-7191-15	0302-7191-15
White	0313-7191-16	0303-7191-16	0312-7191-16	0302-7191-16
Yellow	0313-7191-17	0303-7191-17	0312-7191-17	0302-7191-17

Color-PLUS™ PCL Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCL series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils, mark pens as well as most Inkjet based slide printer, special for Inkjet IPS (Leica) slide printers
Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections




Color-PLUS™ PCL Photosensitive Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCL series microscope slides additionally feature more hydrophobic surface and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, such as pencils, mark pens as well as most laser based slide printer
Recommended Application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4181-02	0303-4181-02	0312-4181-02	0302-4181-02
Green	0313-4181-08	0303-4181-08	0312-4181-08	0302-4181-08
Purple	0313-4181-09	0303-4181-09	0312-4181-09	0302-4181-09
Orange	0313-4181-10	0303-4181-10	0312-4181-10	0302-4181-10
Pink	0313-4181-12	0303-4181-12	0312-4181-12	0302-4181-12
Tan	0313-4181-15	0303-4181-15	0312-4181-15	0302-4181-15
White	0313-4181-16	0303-4181-16	0312-4181-16	0302-4181-16
Yellow	0313-4181-17	0303-4181-17	0312-4181-17	0302-4181-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4191-02	0303-4191-02	0312-4191-02	0302-4191-02
Green	0313-4191-08	0303-4191-08	0312-4191-08	0302-4191-08
Purple	0313-4191-09	0303-4191-09	0312-4191-09	0302-4191-09
Orange	0313-4191-10	0303-4191-10	0312-4191-10	0302-4191-10
Pink	0313-4191-12	0303-4191-12	0312-4191-12	0302-4191-12
Tan	0313-4191-15	0303-4191-15	0312-4191-15	0302-4191-15
White	0313-4191-16	0303-4191-16	0312-4191-16	0302-4191-16
Yellow	0313-4191-17	0303-4191-17	0312-4191-17	0302-4191-17

50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8181-02	0303-8181-02	0312-8181-02	0302-8181-02
Green	0313-8181-08	0303-8181-08	0312-8181-08	0302-8181-08
Purple	0313-8181-09	0303-8181-09	0312-8181-09	0302-8181-09
Orange	0313-8181-10	0303-8181-10	0312-8181-10	0302-8181-10
Pink	0313-8181-12	0303-8181-12	0312-8181-12	0302-8181-12
Tan	0313-8181-15	0303-8181-15	0312-8181-15	0302-8181-15
White	0313-8181-16	0303-8181-16	0312-8181-16	0302-8181-16
Yellow	0313-8181-17	0303-8181-17	0312-8181-17	0302-8181-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8191-02	0303-8191-02	0312-8191-02	0302-8191-02
Green	0313-8191-08	0303-8191-08	0312-8191-08	0302-8191-08
Purple	0313-8191-09	0303-8191-09	0312-8191-09	0302-8191-09
Orange	0313-8191-10	0303-8191-10	0312-8191-10	0302-8191-10
Pink	0313-8191-12	0303-8191-12	0312-8191-12	0302-8191-12
Tan	0313-8191-15	0303-8191-15	0312-8191-15	0302-8191-15
White	0313-8191-16	0303-8191-16	0312-8191-16	0302-8191-16
Yellow	0313-8191-17	0303-8191-17	0312-8191-17	0302-8191-17


Color-PLUS™ PCL Microscope Slides
adaptable to Etching Slide Printers



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The CITOTEST adhesion microscope slides adaptable to etching slide printers always come with a black color coat for contrasting on the opposite side of the marking color-coat side, markable with the conventional labeling system, pencils and mark pens as well as most etching slide printers. The color coating is resistant to most common chemicals and routine stains that are utilized in laboratory.

Recommended application: manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3281-02	0303-3281-02	0312-3281-02	0302-3281-02
Green	0313-3281-08	0303-3281-08	0312-3281-08	0302-3281-08
Purple	0313-3281-09	0303-3281-09	0312-3281-09	0302-3281-09
Orange	0313-3281-10	0303-3281-10	0312-3281-10	0302-3281-10
Pink	0313-3281-12	0303-3281-12	0312-3281-12	0302-3281-12
Tan	0313-3281-15	0303-3281-15	0312-3281-15	0302-3281-15
White	0313-3281-16	0303-3281-16	0312-3281-16	0302-3281-16
Yellow	0313-3281-17	0303-3281-17	0312-3281-17	0302-3281-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3291-02	0303-3291-02	0312-3291-02	0302-3291-02
Green	0313-3291-08	0303-3291-08	0312-3291-08	0302-3291-08
Purple	0313-3291-09	0303-3291-09	0312-3291-09	0302-3291-09
Orange	0313-3291-10	0303-3291-10	0312-3291-10	0302-3291-10
Pink	0313-3291-12	0303-3291-12	0312-3291-12	0302-3291-12
Tan	0313-3291-15	0303-3291-15	0312-3291-15	0302-3291-15
White	0313-3291-16	0303-3291-16	0312-3291-16	0302-3291-16
Yellow	0313-3291-17	0303-3291-17	0312-3291-17	0302-3291-17

Color-PLUS™ PCI Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, pencils, mark pens as well as most thermal transfer based slide printers.

Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (privileged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7161-02	0303-7161-02	0312-7161-02	0302-7161-02
Green	0313-7161-08	0303-7161-08	0312-7161-08	0302-7161-08
Purple	0313-7161-09	0303-7161-09	0312-7161-09	0302-7161-09
Orange	0313-7161-10	0303-7161-10	0312-7161-10	0302-7161-10
Pink	0313-7161-12	0303-7161-12	0312-7161-12	0302-7161-12
Tan	0313-7161-15	0303-7161-15	0312-7161-15	0302-7161-15
White	0313-7161-16	0303-7161-16	0312-7161-16	0302-7161-16
Yellow	0313-7161-17	0303-7161-17	0312-7161-17	0302-7161-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-7171-02	0303-7171-02	0312-7171-02	0302-7171-02
Green	0313-7171-08	0303-7171-08	0312-7171-08	0302-7171-08
Purple	0313-7171-09	0303-7171-09	0312-7171-09	0302-7171-09
Orange	0313-7171-10	0303-7171-10	0312-7171-10	0302-7171-10
Pink	0313-7171-12	0303-7171-12	0312-7171-12	0302-7171-12
Tan	0313-7171-15	0303-7171-15	0312-7171-15	0302-7171-15
White	0313-7171-16	0303-7171-16	0312-7171-16	0302-7171-16
Yellow	0313-7171-17	0303-7171-17	0312-7171-17	0302-7171-17

Color-PLUS™ PCI Microscope Slides



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, pencils, mark pens as well as most inkjet based slide printers.

Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (priviledged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4161-02	0303-4161-02	0312-4161-02	0302-4161-02
Green	0313-4161-08	0303-4161-08	0312-4161-08	0302-4161-08
Purple	0313-4161-09	0303-4161-09	0312-4161-09	0302-4161-09
Orange	0313-4161-10	0303-4161-10	0312-4161-10	0302-4161-10
Pink	0313-4161-12	0303-4161-12	0312-4161-12	0302-4161-12
Tan	0313-4161-15	0303-4161-15	0312-4161-15	0302-4161-15
White	0313-4161-16	0303-4161-16	0312-4161-16	0302-4161-16
Yellow	0313-4161-17	0303-4161-17	0312-4161-17	0302-4161-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-4171-02	0303-4171-02	0312-4171-02	0302-4171-02
Green	0313-4171-08	0303-4171-08	0312-4171-08	0302-4171-08
Purple	0313-4171-09	0303-4171-09	0312-4171-09	0302-4171-09
Orange	0313-4171-10	0303-4171-10	0312-4171-10	0302-4171-10
Pink	0313-4171-12	0303-4171-12	0312-4171-12	0302-4171-12
Tan	0313-4171-15	0303-4171-15	0312-4171-15	0302-4171-15
White	0313-4171-16	0303-4171-16	0312-4171-16	0302-4171-16
Yellow	0313-4171-17	0303-4171-17	0312-4171-17	0302-4171-17

Color-PLUS™ PCI Photosensitive Microscope Slides




Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.

Markable with the conventional labeling system, pencils, mark pens as well as most laser based slide printers.

Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (priviledged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8161-02	0303-8161-02	0312-8161-02	0302-8161-02
Green	0313-8161-08	0303-8161-08	0312-8161-08	0302-8161-08
Purple	0313-8161-09	0303-8161-09	0312-8161-09	0302-8161-09
Orange	0313-8161-10	0303-8161-10	0312-8161-10	0302-8161-10
Pink	0313-8161-12	0303-8161-12	0312-8161-12	0302-8161-12
Tan	0313-8161-15	0303-8161-15	0312-8161-15	0302-8161-15
White	0313-8161-16	0303-8161-16	0312-8161-16	0302-8161-16
Yellow	0313-8161-17	0303-8161-17	0312-8161-17	0302-8161-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-8171-02	0303-8171-02	0312-8171-02	0302-8171-02
Green	0313-8171-08	0303-8171-08	0312-8171-08	0302-8171-08
Purple	0313-8171-09	0303-8171-09	0312-8171-09	0302-8171-09
Orange	0313-8171-10	0303-8171-10	0312-8171-10	0302-8171-10
Pink	0313-8171-12	0303-8171-12	0312-8171-12	0302-8171-12
Tan	0313-8171-15	0303-8171-15	0312-8171-15	0302-8171-15
White	0313-8171-16	0303-8171-16	0312-8171-16	0302-8171-16
Yellow	0313-8171-17	0303-8171-17	0312-8171-17	0302-8171-17


Color-PLUS™ PCI Microscope Slides
adaptable to Etching Slide Printers



Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

The Color-PLUS™ PCI series microscope slides are especially designed for the automatic IHC staining including the Roche Ventana automated IHC Stainers. The microscope slides feature extraordinary adhesion ability and more hydrophilic surface, and always come with a silky color coating at one surface End that is resistant to all common chemicals and routine stains that are utilized in laboratory.


Markable with the conventional labeling system, pencils, mark pens as well as most etching based slide printers.
Recommended Application: automatic IHC staining with the Roche Ventana automated IHC Stainer (priviledged use) , manual IHC staining, automatic IHC staining with the Leica and Dako automated IHC Stainers and routine H&E staining of detachable tissue sections.



Liquid-Base Cytology Microscope Slides

Dimensions: 25.0mm×75.0mm Thickness: 1.0-1.2mm

Supplied for liquid-base cytology, e.g., TCT & LCT slide preparation.



- Ground Edges, 45°clipped corners
- Printed Circle or Bracket optional
- Coated and non-coated available
- Frosted area of approx. 20mm depth, displays at one surface end

Packing:
50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton
72pcs/box, 1,440pcs/carton

50/PK Packing: 50pcs/box, 1,000pcs/promo-pk, 2,000pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3261-02	0303-3261-02	0312-3261-02	0302-3261-02
Green	0313-3261-08	0303-3261-08	0312-3261-08	0302-3261-08
Purple	0313-3261-09	0303-3261-09	0312-3261-09	0302-3261-09
Orange	0313-3261-10	0303-3261-10	0312-3261-10	0302-3261-10
Pink	0313-3261-12	0303-3261-12	0312-3261-12	0302-3261-12
Tan	0313-3261-15	0303-3261-15	0312-3261-15	0302-3261-15
White	0313-3261-16	0303-3261-16	0312-3261-16	0302-3261-16
Yellow	0313-3261-17	0303-3261-17	0312-3261-17	0302-3261-17

72/PK Packing: 72pcs/box, 1,440pcs/carton



MARKING AREA	GROUND EDGES, 90°CORNERS		GROUND EDGES, 45°CORNERS	
	CSW	SUPER GRADE	CSW	SUPER GRADE
Blue	0313-3271-02	0303-3271-02	0312-3271-02	0302-3271-02
Green	0313-3271-08	0303-3271-08	0312-3271-08	0302-3271-08
Purple	0313-3271-09	0303-3271-09	0312-3271-09	0302-3271-09
Orange	0313-3271-10	0303-3271-10	0312-3271-10	0302-3271-10
Pink	0313-3271-12	0303-3271-12	0312-3271-12	0302-3271-12
Tan	0313-3271-15	0303-3271-15	0312-3271-15	0302-3271-15
White	0313-3271-16	0303-3271-16	0312-3271-16	0302-3271-16
Yellow	0313-3271-17	0303-3271-17	0312-3271-17	0302-3271-17

Microscope Slides with Printed Circles

Compatible with single or double Cellfunnels for use in cyto centrifuges. Recommended for Thermo Scientific™ Shandon™ EZ Cytotunnel™ or other brands cytofunnels which have the same design.

	MARKING AREA	DESCRIPTION	50PCS		72PCS	
			CSW	SUPER GRADE	CSW	SUPER GRADE
A	White	Single printed circle Dia.: 13mm	0370-0003	0370-0103	0370-0203-16	0370-0303-16
	White	Single printed circle Dia.: 13mm, coated	0370-0001	0370-0101	0370-0201-16	0370-0301-16
B	White	Double printed circles Dia.: 13mm	0370-0004	0370-0104	0370-0204-16	0370-0304-16
	White	Double printed circles Dia.: 13mm, coated	0370-0002	0370-0102	0370-0202-16	0370-0302-16
C	White	Double printed circles Dia.: 13mm	0370-0020	0370-0120	0370-0220-16	0370-0320-16
	White	Double printed circles Dia.: 13mm, coated	0370-0021	0370-0121	0370-0221-16	0370-0321-16

Microscope Slides with Printed Bracket

Specifically designed for use with the ThinPrep® and other brand processors.

	MARKING AREA	DESCRIPTION	50PCS		72PCS	
			CSW	SUPER GRADE	CSW	SUPER GRADE
D	White	Printed Bracket Dia.: 36mm, coated	0370-0013-16	0370-0113-16	0370-0213-16	0370-0313-16

