

## Blood Collection Minitor

### BMC-1200A

- ◆Function of load-bearing, Blood bag weight deducted and Liquid production control.
- ◆Functions of pallet swing, alarm, proximity control, swing stop and swing automatic recovery
- ◆Eliminates the man-made error and provides convenience for the production of plastic bags
- ◆Suitable for various blood sampling mechanism

### BMC-1200B

- ◆For blood stations, hospitals and other units to separate the use of plasma
- ◆After bags in the blood is layered
- ◆Points pulp clamp pressure can make the different composition of blood top-down being squeezed out, heating constant temperature function



BMC-1200A



BMC-1200B

Model	BMC-1200A	BMC-1200B
Range	0-1200ml	0-1200ml
Division value	1ml	1ml
Liquid production rate	0.5--3ml/s	0.5--3ml/s
Liquid production error	± 5%	± 5%
Liquid specific gravity	1.05g/mL	1.05g/mL
Oscillation frequency	30-32rpm	30-32rpm
Swing angle	13±2	13±2
Power	40W	40W
Alarm and block	Sound light alarm and block	Sound light alarm and block
Block applied piping	--	external diameter: Ø3~Ø6mm, thickness: 0.5mm
Block thrust	--	≤14N
Voltage	110V 60Hz/ 220V 50Hz	110V 60Hz/ 220V 50Hz
Dimension(mm)	220x170x160	275x230x210
Package (mm)	290x240x210	350x300x260
G.W/N.W.	3/2.5kg	5.5/5kg