



Fetal Doppler

8 C/8 D

8 series ultrasonic fetal doppler meets the fetus daily check and routine examination at home, clinic, community and hospital.



PC Software

Features

- High-fidelity, crystal clear sound
- Earphone and speaker are feasible
- Display with backlight or color LCD
- USB functions to transmit data
- Multiple display modes
- Both audio and light alarm
- Multiple work modes
- Special antiskid function



Fetal Doppler




Technical Specification

- **Ultrasound frequency:** 2MHz/2.5MHz/3MHz
- **Ultrasound intensity:** <10mW/cm²
- **Power supply:**
DC Ni-Mh rechargeable battery
AC 220/110V, 50/60Hz
(for D and C)
- **Display:** 65mm×50mm LCD
- **FHR resolution:** 1bpm
- **FHR accuracy:** ±2bpm
- **Power consumption:** <12W
- **FHR Measuring range:** 50~230bpm
- **Dimension:** 220mm×220mm×140mm
- **Weight:** 2.1Kg

PC software

Data transmit, edit, store, analyse, print, etc.



Standard Accessories

2MHz Probe



Adapter



Gel



Optional

2.5/3MHz Probe



PC Software



Carry Case



Trolley



Print function



All specifications are subject to change without prior notice.